"Exciting and detecting waves in micro and nanoscale systems by Picosecond Ultrasonics"





Jean-François ROBILLARD jrobilla@email.arizona.edu University of Arizona - ERC TeleSeminar August 12th 2010

Who am I?

Master's Physics PhD PhD Viversité Lille1 Vices et remove Picosecond Ultrasonics Today's talk Postdoc Magneto-elasticity

> Research Associate Materials Science & THE UNIVERSITY OF ARIZONA. Engineering (Pr. Deymier) NanoPhononic crystals





jrobilla@email.arizona.edu



What am I interested in ?

•Sound waves (Guiding, Use for metrology, Ultra-high frequencies)

• Phononic Crystals (Design, Functionalities, Tunability)

•Nanoscale materials (Thermal and acoustic properties)

Couplings (Magneto-acoustics)

Picosecond Ultrasonics





Picosecond Ultrasonics

- Optical pump-probe experiment
- Uses femtosecond laser pulses
- To generate/detect acoustic waves

Scales the sonar principle to the µm/nm range *i.e.* GHz/THz frequency range





Setup – Main parameters





Laser wavelength λ Second Harmonic Generation RED Pump BLUE Probe





Another example







Orders of magnitude

Acoustic pulse width = 2 * light absorbtion depth

2*10= 20 nm in metals

Duration of pulse = pulse width/velocity

20 nm/(5 nm/ps) = 4 ps

• Frequency content of acoustic pulse = 1/duration

1 / 4ps = 250 GHz

• Max thickness of sample = $\frac{1}{2}$ velocity * time between two light pulses $\frac{1}{2} * 5(nm/ps) * 13000 \text{ ps} \approx 30 \mu \text{m}$





A wealth of detection mechanisms...



Attenuation in SiO2 @ 250 GHz

Thickness $@\pm 2 \text{ nm}$



A. Devos et al. App. Phys. Lett., 86, 211903, (2005)A. Devos et C. Lerouge. Phys. Rev. Lett., 86, 2669, (2001)

A. Devos, J.-F. Robillard et al. Phys. Rev. B, 74, 064114, (2006)



Summary - Abilities and Uses of Picosecond Ultrasonics

10

Formation Transfert

•	Measures time	W
	Echos	Al
	Brillouin oscillations period	Cr Mo
•	To get	Ni
	Thickness (30 Å – 10 µm)	Ta
	Acoustic properties (Sound velocity, Elastic modulus, Attenuation)	Au
•	Solid state physics	TiNi Si
	Semi-conductors physics, Attenuation in glasses	Ge SiGe
	Nanotechnologies, Superlattices, quantum dots, colloïds	GaAs GaP
•	Thin films metrology	GaN AlGaN
	Sound velocity, thickness, adherence, multilayers stacks	SiO2
	Wide variety of materials, sound waves propagate in every solid !	SiOC
	(Metals, SC, Dielectrics, Piezo)	AIN PZT
	Industry use PU in-line (Metapulse Rudolph techs)	SrTiO3



Exciting and detecting waves in lattices of nanoscale objects with Picosecond Ultrasonics







Why?

- Nano-objects vibrations
 - GHz THz
 - Few studies
- Hypersonic crystals
 - Acoustic devices
 - Frequency ↔ lattice constant
- Multi-«onic» Crystals
 - Phonons
 - Photons
 - Plasmons



AI

$$f = \frac{c}{d}$$
 f = 32 GHz
T = 30 ps
200 nm
 $f = 32 \text{ kHz}$
T = 30 µs

$$\mu m \rightarrow GHz$$



PU and nano-objects

PU is suitable for the study of eigenmodes



C. Voisin, D. Christofilos, N. DelFatti and F. Vallée Physica B **316–317**, 89–94 (2002)



H.-N. Lin, H. J. Maris et al. J. Appl. Phys. **73**, 37-45 (1993)





What about periodic, phononic structures



How ? E-beam lithography

- Electrons print photo-resist
- Electrons wavelength $\lambda_e << 1$ nm
- High resolution ~ a few nm

LEICA EBPG5000 plus

Voltage 20 – 100kV Beam Current 100pA – 200nA Spot diameter < 10nm









TUCSON ÁRIZONA







Image: Source A = InLens Signal A = InLens Signal = 0.3846 WD = 2 mm Image: Source K x Signal B = SE2 Mixing = Off WD = 2 mm

Au / Ti a = 600 nm d = 200 nm Substrate Si



Al a = 700 nm d = 200 nm Sur Al h = 100 nm Substrat pyrex THE UNIVERSITY OF RIZONIA

TUCSON ARIZONA

Samples



Al a = 300 nm d = 100 nm Sur Al h = 100 nm Substrate fused silica



Al a = 700 nm d = 200 nm Sur Al h = 100 nm Substrat pyrex



Al a = 400 nm d = 200 nm Substrat eSi



Pt a = 1000 nm d = 500 nm Sur Al h = 100 nm Substrat pyrex

Individual Modes





Individual modes – Experimental evidence

The probe collects signals from a large number of cubes...



Size d \propto lattice constant \rightarrow elastics Magnitude \rightarrow sound velocity

Elastic Contribution

J.-F. ROBILLARD, A. DEVOS and I. ROCH-JEUNE Physical Review B 76, 092301 (2007)



THE UNIVERSITY OF ARIZONA TUCSON ARIZONA

Individual modes – Experimental evidence



- Period do not depends on a
- Signal ∞ number of cubes per unit surface



Signal comes from cubes !



Modeling individual modes



- Complex objects
- Strong couplings (high attenuation)



Collective Modes





Collective modes



Same sample, two probe wavelengths A new signal •Lower frequency •Higher magnitude •Longer lifetime

TUCSON ARIZONA

Detection criterion:





ARIZONA TUCSON ARIZONA

Formation Transfert

Method





Transfert

Modeling

 $\omega(\vec{k}) = c(a) \times \left\| \vec{k} \right\|$

Mode velocity -Curvature -Limit velocity

Wavevectors -Branches -Factors





Modèle élastique – Vecteurs d'onde

• The generated mode is periodic $\forall \vec{A}$ Direct lattice $\eta_i(\vec{r},t) = \eta_i(\vec{r}+\vec{A},t)$

$$\eta_i(\vec{r},t) = \eta_{0i} \exp(i\vec{k}.r - \omega t)$$
$$\vec{k}.\vec{A} = 2n\pi$$

Wavevectors belong to the reciprocal lattice Exemple : Square lattice

$$\|\vec{k}_{10}\| = \frac{2\pi}{a} \qquad \|\vec{k}_{20}\| = \frac{2\pi}{a} \times 2$$
$$\|\vec{k}_{11}\| = \frac{2\pi}{a} \times \sqrt{2} \qquad \|\vec{k}_{21}\| = \frac{2\pi}{a} \times \sqrt{5}$$
Each wavevector gives a branch

In-plane propagation, several directions





Modeling – Mode velocity

Lattice cell = harmonic oscillator

$$\omega = \sqrt{\frac{K}{M}} = \sqrt{\frac{K}{\rho(a^2h + d^3)}} = c(a) \|\vec{k}_{10}\|$$
$$c(a) = \frac{1}{2\pi} \sqrt{\frac{K}{\rho}} \sqrt{\frac{a^2}{a^2h + d^3}}$$

Limit speed = « without cubes »

$$\lim_{a\to\infty} c(a) = c_0$$

$$c(a) = c_0 \sqrt{\frac{a^2 h}{a^2 h + d^3}}$$





Modeling - Discussion







Modeling - Discussion







Model vs experiments

Lattice symetry





THE UNIVERSITY OF ARIZONA TUCSON ARIZONA



Model vs experiments



Size of cubes d



Layer thickness h



Layer / substrate system





Use for thin films metrology





Use for metrology

- PU « as usual »
 - Longitudinal waves
 - Thickness, longitudinal speed and density
- «Nano Periodic transducer»
 - Surface waves
 - Surface waves velocity
- Advantages
 - Complete characterization
 - All materials
 - Standard PU setup
 - Thin films < 100 nm



Use for metrology

- Test sample SiO₂
- Two sound velocities
 - Longitudinale c_l
 - Rayleigh c_{Rayleigh}
- Deduce:

Young modulus E = 72 ± 3 GPa Poisson ratio $\upsilon = 0.16 \pm 0.08$

Agreement with litterature





P. A. MANTE, **J. F. ROBILLARD** and A. DEVOS, Applied Physics Letters **93**, 071909 (2008) 37

Transfert

Conclusions

Picosecond Ultrasonics

- A nano to micro GHz to THz sonar
- Useful for solid state physics and thin films metrology

Lattices of nanoscale objects

- Individual modes
 - Complex objects (many modes, strong couplings)
- Collective modes
 - Technology, e-beam
 - Efficient analytic model
 - Opens up another metrologic capacity

Perspectives

- Semiconductors lattices, namowires, nanocrystals
- Anisotropic materials
- Thermal properties
- Higher frequency phononic crystals... « Thermonic crystals »

Thanks for your attention !



jrobilla@email.arizona.edu



Acknowledgements



• A. Devos, I. Roch R. Côte, P. Emery and P.A. Mante



 \bullet

B. Perrin, T. Bienville, L. Beliard

ullet



THE UNIVERSITY

P. A. Deymier, K. Muralidharan, W. Beck, S. Seraphin, D. Barker, J. Bucay and N. Swinteck

THE UNIVERSITY OF ARIZONA TUCSON ARIZONA US National Science Foundation International Lab MATEO (UofA / CNRS)

